

WHAT IS IT?

RIGENERA is a device specially developed for the clinical application of regenerative procedures using autologous micrografts. In the same session, the patient is a donor and receiver of micrografts, which allows the receptor area to benefit from the regenerative activity of the cellular content and the growth factors extracted from the donor area.

COMPONENTS

RIGENERA is a patented technology for the safe preparation of micrografts through mechanical disintegration and automatic filtering of any solid tissue.



RIGENERA MACHINE

Is a compact unit that works with the disposable Rigeneracons and operates at a constant speed of approximately 80 rpm.




RIGENERACON

Is a disposable medical device with CE and FDA certification. It contains a stationary stainless steel grid of 2.5ml capacity and it is made up of 100 hexagonal holes. Around each hole there are 6 microblades, with a total of 600 microblades designed to efficiently cut hard and soft tissues without damaging their structure.

WHAT IS IT BASED ON?


It is based on the Tehory of **MICROGRAFT TISSUE**.
By using **HEALTHY HOMOLOGOUS TISSUE (AUTOLOGOUS)** + **RIGENERA** the result is **REGENERATION OF DAMAGED TISSUE**.

CONVENTIONAL MICROGRAFT



1:1 PROPORTION
With 1 cm² of healthy tissue (biologically homologous) it regenerate 1 cm² of damaged tissue.

NEW THEORY OF MICROGRAFT




1:20 PROPORTION
With 1 cm² of healthy tissue (biologically homologous) it regenerate 20 cm²

IF WE USE HOMOLGY REGENERATION

When we extract **healthy connective tissue HOMOLOGOUS** to the damaged tissue that we want to regenerate, **we take advantage of 95% of its:**

-  DIFFERENTIATION CAPACITY
-  POTENTIALITY
-  EFFICIENCY


***HOMOLOGOUS TISSUE**
■ REGENERATION
■ REPAIR



***Homologous Tissue:** Connective tissue biologically equal to the one that needs to be regenerated.

THE RESULT IS 80% REGENERATION AND 20% REPAIR

REGENERATING IS DIFFERENT THAN REPAIRING



Skin Rejuvenation

Micrografting helps the dermatologist, aesthetic physician or plastic surgeon to obtain better results in its multiple procedures for the treatment of elastosis, skin aging and healing, as well as a much faster recovery period for the patient.

WHAT DOES THE MICROGRAFTING OFFER?


- Contributes to the treatment carried out, growth factors as well as a tissue regeneration in the face of the controlled aggression of the done treatment.
- It is indicated in what degree of skin aging as an adjuvant to the treatment chosen by the specialist.
- It also promotes regeneration and controls the periods of inflammation after performing a surgical suture.
- It helps all treatments to improve:



Skin texture



Reduce depth wrinkles



Decrease sunspots and aging



ADVANTAGES

Without cellular cultivation

Without chemical and biological treatment

No side effects (autologous)

Fast and simple procedure (30 minutes)

No hospitalization / operating room needed

Only a local anesthetic is necessary

No extra specialized equipment needed



PROCEDURE

Requirement: Once the treatment chosen for skin rejuvenation has been carried out, the micrografting procedure is implemented.

1

Perform chosen procedure for skin rejuvenation

De-epidermalization with laser, dermabrasion etc

2

The donor sample is obtained by means of a 4 mm punch of the mastoid hairy area

With 2 punch of 4 mm is enough to treat an area of 20 cm²

3

Process the micrograft by RIGENERA MACHINE

Get 1.2 ml. of solution to inject after 2 minutes of processing the biopsied sample.

4

Inject the solution intradermally through mesotherapy, creating papillae in the treated area.

Use 1 ml syringe luer lock and 30 G needle.

BIBLIOGRAPHY

- Letizia Trovato, Manuela Monti, Claudia del Fante, Marila Cervio, Mila Lampinen, Lucia Ambrosio, Carlo Alberto Redi, Cesare Perotti, Esko Kankuri, Gennaro Ambrosio, Ruggero Rodriguez Y.Baena, Giuseppe Pirozzi and Antonio Graziano. *A New Medical Device Rigeneracons Allows to Obtain Viable Micro-Grafts From Mechanical Disaggregation of Human Tissues*. J Cell Physiol 2015;230:2299–303.
- Marco Marcarelli, Letizia Trovato, Elvio Novarese, Michele Riccio and Antonio Graziano. *Rigenera protocol in the treatment of surgical wound dehiscence*. 2016. International Wound Journal ISSN 1742-4801
- Fabiano Svolacchia, Francesco De Francesco, Letizia Trovato, Antonio Graziano, Giuseppe A. Ferraro. *An innovative regenerative treatment of scars with dermal micrografts*. J Cosmet Dermatol 2016. DOI: 10.1111/jocd.12212
- Elisabetta Baglioni, Letizia Trovato, Marco Marcarelli, Ambra Frenello and Maria Alessandra Bocchiotti. *Treatment of Oncological post-Surgical Wound Dehiscence with Autologous Skin Micrografts*. Anticancer Research 36: 975-980 (2016)
- Valeria Purpura, Elena Bondioli, Antonio Graziano, Letizia Trovato, Davide Melandri, Martina Ghetti, Andrea Marchesini, Maria Gabriella Cusella De Angelis, Laura Benedetti, Gabriele Ceccarelli, Michele Riccio. *Tissue Characterization after a New Disaggregation Method for Skin Micro-Grafts Generation*. (109), e53579,doi:10.3791/53579 (2016).